

Product Name :

Sterilising Fluid

Product Code :

SLE/SAH/30

**Description :**

An essential sterilising fluid for protective eyewear and equipment.

Kills 99.9% of germs including bacteria, fungi and viruses.

Soak equipment for 15 minutes in solution with 30 ml (one capful) of fluid in 5 litres of water, or a dilution of 0.6%.

Effective on Rotavirus, Coronavirus, MRSA, swine flu and more.

Sterilising solution will be effective for up to 24 hours.

Leave equipment to dry, no need to rinse.

Available in 500ml, 1000ml and 5 L bottles.

Contact Science Lab for your Educational School Science Lab Equipments. We are best mathematics lab equipments suppliers, mathematics lab kit, mathematics lab kits india, maths kits ,maths lab equipment manufacturers in india, maths lab equipment suppliers.