



THE LAB EQUIPMENT

Email : [sales@thelabequipment.com](mailto:sales@thelabequipment.com)

Phone: +91-8759245890

**Product Name :**

15-Lead ECG Placement Trainer

**Product Code :**

SLE/A/141



## **Description :**

Teaches up to 15-lead ECG electrode placement anatomically and the accuracy of the electrode placement.

Students learn the placement of the electrodes on the adult trainer using anatomical landmarks, such as intercostal spaces, midclavicular line, anterior axillary line, midaxillary line, and scapula. The trainer does not provide ECG output signals, but simulations of rhythms and hands-free defibrillation can be performed by adding any ECG arrhythmia simulator directly to an ECG monitor.

Connection sites for 4 limb leads and V1 through V9, with the ability to attach right- or left-sided electrodes.

**Website:** <https://thelabequipment.com>, **Email:** [sales@thelabequipment.com](mailto:sales@thelabequipment.com)

plot. no.296. Sector 5. Vaishali. Ghaziabad. Uttar Pradesh 201010